DE	Rückenschmerzgürtel Gebrauchsanweisung	2
EN	Back pain belt Instructions for use	12
FR	Ceinture lombaire Mode d'emploi	21
ES	Cinturón terapéutico lumbar Manual de instrucciones	31
IT	Cintura per dolori lombari Istruzioni per l'uso	41
TR	Sırt ağrısı kemeri Kullanma Talimatı	51
RU	Электростимулятор для спины Инструкция по применению	60
PL	Pas likwidujący ból pleców Instrukcja obsługi	70

Inhalt

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Germany

Contents

1. About the unit.....12	7. Cleaning and care of the unit and belt18
2. Signs and symbols.....12	8. Disposal.....18
3. Important notes.....13	9. Troubleshooting.....19
4. Unit description.....16	10. Technical specifications19
5. Using the unit for the first time and changing the batteries.....16	11. Notes on electromagnetic compatibility.....20
6. Operation.....17	12. Warranty/service.....20

Included in delivery

- Back belt with control unit EM 38
- Extension belt
- Batteries, 3 x 1.5 V AAA (LR03, micro)
- This instruction manual

1. About the unit

Dear valued customer

Thank you for choosing one of our products.

Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

With kind regards
Your Beurer team

Pain relief on lower back with the EM 38 stimulation current unit

How does stimulation current work?

The back belt works on the basis of electrical nerve stimulation (TENS).

The more precise term is transcutaneous (meaning “through the skin”) electrical nerve stimulation (TENS). TENS is approved as a clinically proven, effective, medication-free method of treating certain types of pain which, when used properly, is also free of side effects. It is also approved as a means of simple self-treatment.

The pain-relieving or pain-suppressing effect is caused by suppressing the transmission of pain in the nerve fibres (particularly through high-frequency impulses) and increasing the release of the body’s own endorphins, which decrease the sensation of pain by their effect in the central nervous system. This method is clinically tested and approved.

Every condition that calls for the use of TENS must be clarified by your attending physician. He or she will also give you information about the specific benefits of a TENS self-treatment.

2. Signs and symbols

The following symbols appear in these instructions.

	Warning	Warning instruction indicating a risk of injury or damage to health.
	Important	Safety note indicating possible damage to the unit/accessories.
	Note	Important information to be noted.

The following symbols are used on the identification plate.

	Type BF applied part
	Refer to instructions for use
	Devices that can emit effective output values above 10mA or 10V, averaged over every 5 second interval.
	Manufacturer
	Keep dry
	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.
	Certification symbol for products that are exported to the Russian Federation and CIS countries.

3. Important notes

Safety notes

Warning

- Use the belt only:
 - on humans,
 - for external use,
 - for the purpose for which it was developed and in the manner specified in these instructions for use.
- Any improper use can be dangerous!
- In acute emergencies, first aid takes priority.
- This unit is not intended for commercial or clinical use, but only for individual, private household use!
- Before use, ensure that there is no visible damage to the unit or accessories. When in doubt, do not use the unit and contact your dealer or the customer service address provided.
- If you have health concerns of any kind, contact your general practitioner!
- Use the unit in the lower back area only! Using it on other body parts can cause severe health problems.
- A slight reddening of the skin after use is normal and disappears after a short time. Refrain from using the unit until the reddening has gone away.
- If skin irritation results from a long therapy time, select a shorter application time.
- In case of more severe skin irritation, discontinue treatment and consult a physician.
- This unit is not intended for use by children or persons with physical, sensory (e.g. insensitivity to pain) or mental impairments or those who do not know how to use it due to a lack of experience or knowledge, unless they are supervised by a person responsible for their safety or have received instructions from them on how to use the unit.
- Keep children away from packaging materials (risk of suffocation).
- Do not use attachments that are not recommended by the manufacturer.

Safety precautions

Warning

- Carry out the first few minutes of treatment while sitting or lying down, so that if you experience a rare vasovagal reaction (i.e. you feel faint), you are not at unnecessary risk of injury. If you start to feel faint, shut off the unit immediately and elevate your legs (for approx. 5 – 10 min.).
- The treatment should be pleasant. If the unit does not function correctly or if you begin to feel unwell or feel

- pain, stop using the unit immediately.
- Take off the belt only when the unit is switched off!

Notes about the electrodes

Important

- Never apply the electrodes to injured skin areas.
- Max. recommended initial value for electrodes: 5 mA/cm².
- Effective current density levels over 2 mA/cm² require increased attention on the part of the user.

Warning

To prevent harm to your health, we urgently recommend that you do not use the unit in the following situations:

- Do not use the unit if you use a pacemaker or other implants, such as an insulin pump or metallic implants.
- If you have a high fever (e.g. > 39 °C).
- If you have known or acute arrhythmias or other problems involving cardiac stimuli and their transmission to the heart.
- If you have seizure disorders (such as epilepsy).
- If you are pregnant.
- If you have cancer.
- After operations in which increased muscle contractions could interfere with the healing process.
- On acutely or chronically diseased (injured or inflamed) skin, e.g. painless or painful inflammation or reddening.
- If you have rashes (e.g. allergies), burns, bruises, swelling and open wounds or those in the process of healing.
- On surgical scars in the process of healing.
- When simultaneously connected to a high-frequency surgical instrument. In this case, burns under the stimulation current fields can result.
- When under the influence of pain medication, alcohol or sleep aids.
- Whenever doing anything in which an unforeseen reaction (such as increased muscle contraction despite low intensity) can be hazardous, for example while driving a car or operating machinery.
- On sleeping persons.
- For acute or chronic diseases of the gastrointestinal tract.
- Do not use this unit simultaneously with other units that transmit electrical pulses to your body.
- The unit is suitable for private, individual use.
- For hygienic reasons, the belt may be used by one person only.
- Ensure that during stimulation, no metallic objects can contact the electrodes, as otherwise localised burns can occur.

The unit must not be used in the following situations:

- In the head area: Here, it can trigger cramping.
- In the neck/carotid artery area: Here, it can trigger cardiac arrest.
- In the throat and larynx area: Here, it can trigger muscle cramps that can cause suffocation.
- In the ribcage area: Here, it can increase the risk of ventricular fibrillation and cause cardiac arrest.

Safety precautions

- Do not use the unit in close proximity to shortwave or microwave units (e.g. cell phones), as these can cause fluctuations of the initial values of the unit.
- Do not use while sleeping, driving a motor vehicle or simultaneously operating machinery.
- Never immerse the unit in water or other liquids.
- Do not use in close proximity to highly inflammable materials, gases or explosives.

Consult a doctor before use in the following situations:

- Acute diseases, particularly in case of suspected or known blood coagulation disorders, susceptibility to thromboembolitic disease or malignant growths.
- If you have diabetes or other diseases.
- Chronic pain conditions that have not been clarified, regardless of the region of the body.
- Any sensitivity disorders with reduced sensation of pain (such as metabolic disorders).
- Simultaneously carried out medical treatments.
- Discomfort experienced as a result of the stimulation treatment.
- Persistent skin irritation under the electrodes.

Warning

Using the unit is not a substitute for consultation or treatment from a physician. Whenever you have any pain or illness, always contact your physician first.

Before using the unit for the first time

Important

- Before you use the unit for the first time, remove all packaging materials.
- Switch the unit off immediately if it is defective or malfunctioning.
- Never wear the belt with bare metal electrodes. Using the belt without electrode covers can cause injuries.

Important

- The manufacturer shall not be held liable for damage or injuries caused by improper or incorrect use.
- Protect the unit from dust, dirt and moisture.
- If the unit has fallen or been dropped, exposed to extreme moisture or incurred other types of damage, it must no longer be used. Do not expose the unit to high temperatures or direct sunlight.
- Never open or attempt to repair the unit by yourself, as otherwise proper function is no longer guaranteed. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.


Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
-  **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

4. Unit description

Overview of the control unit

Item	Item Designation
1	On/off/pause button
2	Intensity adjustment ▲ increase ▼ decrease
3	Program/locking key
4	Battery compartment
5	Display with LCD display


Item	Item Designation
6	Battery status, appears when batteries are used up
7	Intensity, level 0 – 20
8	Pause symbol, flashes when activated
9	TIMER: Remaining time of the active programs in minutes
10	Keylock active
11	Active program

5. Using the unit for the first time and changing the batteries

Putting on the belt

Before each use, ensure that your skin is clean and free of creams or oils. Only in this way the unit can work optimally.

- Lay out the belt before you with the electrodes facing you.
- Moisten the electrodes and your back with water. Ensure that the electrode surfaces are moistened evenly to avoid localised voltage peaks when they are applied to your skin.
- Fit the belt around your waist so that the vertical seams **A** are at the same place on each side of your body. The Beurer logo **B** is then in the middle of the spinal column.

i See the diagram.

When positioned correctly, the electrodes **C** are slightly to the left and right of your spinal column.

- Then, close the Velcro fastener.

i Ensure that the belt fits tightly enough to provide a good contact between the electrodes and your skin, but does not constrict.


i **Extending the belt**
If the belt is not adequate for your waistline, you can use the extension belt provided. Depending on your waistline, the holder for the control unit may be away from the centre of your body. It is important that the belt is positioned correctly, therefore the electrodes **C** must be to the left and right of the spinal column.

Inserting the control unit

- Insert the control unit into the holder on the belt.
- The operating keys should be facing outwards and the display should be facing upwards and towards you.
- If inserted correctly, the control unit clicks into the holder.

Automatic shutoff



If the unit is not used after being switched on, it switches off automatically after 3 minutes.

i Changing the batteries

If the battery status indicator [6] flashes, it is time to change the batteries. To insert or replace batteries, you have to remove the control unit from the holder of the belt.

- Release the magnetic connection by gently pulling the control unit away from the belt. You can feel the magnetic connection of the two contact buttons unsnap.
- Pull the unit upwards and out of the holder.
- Open the battery compartment [4] by pushing the cover downwards at the arrow mark.
- Insert 3 x 1.5 V batteries, type AAA (LR03 Micro). Ensure the correct polarity. See the diagram inside the battery compartment.
- Do not use any rechargeable batteries!

6. Operation

General information about use

The back belt is intended for individual treatment of pain in the lower back area.

i Do not switch on the unit until the belt is fitted correctly. See item 5.

Switching on the unit

Press and hold down the on/off button [1] until a brief signal sounds and the LCD display [5] switches on. When the unit is started for the first time, Program A is activated automatically.

Selecting the program

Press the program button P [3] to switch through the programs. The following programs are available:

Program	Frequency	Time
A	4 Hz –110 Hz (3 phases)	30 min.
B	4 Hz	25 min.
C	2 Hz (burst)	25 min.
D	100 Hz	25 min.

i **Note:**
For program A, you can feel the effect become stronger when switching from phase 1 to phase 2 (after approx. 10 minutes). This is normal and intentional. If the intensity seems too high, decrease it using the ▼ intensity button [2].

i Note:

If you change the program during the stimulation (for example from A to B), the initial intensity in the new program increases incrementally up to the previously set intensity. This can be stopped if necessary by pressing the ▼ intensity button [2] for two seconds or by switching off the unit by pressing the on/off button [1] for two seconds.

Stopping impulses that are too strong

You can decrease the intensity at any time or switch off the unit by pressing the on/off button [1] (approx. 2 seconds).

Adjusting the intensity

Press the ▲ intensity button [2] to increase the intensity incrementally or press the ▼ intensity button [2] to lower the intensity. You can adjust the intensity in 20 Increments. Depending on the intensity level, you feel a prickling at first that can escalate to a muscle contraction.

i

Select a setting that is within your comfort level.

Preventing unwanted pulse changes

To prevent an accidental increase in intensity during a treatment session, simply activate the keylock. To do so, press the program selection button P [3] for approx. 2 seconds. An acoustic signal sounds and the “” symbol appears in the display[5].

To cancel the keylock, press and hold down the P [3] button for approx. 2 seconds.

Reacting to unpleasant sensations

If you feel a stinging or itching sensation, switch off the unit and proceed as follows:

- Check whether the electrodes are working properly or are defective.
- Check whether the round cover is still present over the electrode connections.
- Take off the belt and moisten it completely once again.
- When refitting the belt, ensure good skin contact and good moistening.

Inserting a pause

If you want to take a break or pause a program, press the pause button [1] briefly. An acoustic signal sounds and the “**II**” pause symbol flashes in the display [5]. You can resume the program by pressing the button [1] again.

Contact detection

If the electrodes have no contact with the body, the intensity is set automatically to zero. This prevents unwanted stimulation. The intensity cannot be increased without contact with the body.

Memory function

The unit stores the most recently configured program.

After a battery change, the unit starts with the first program again.

7. Cleaning and care of the unit and belt

Cleaning the unit

- Remove the batteries from the unit before each cleaning.
- Clean the unit with a soft cloth moistened slightly with water after use. If the dirt is heavier, you can also moisten the cloth with a mild soap solution.
- Ensure that no water penetrates the unit. If this should occur, do not use the unit until it is completely dry.
- Do not use any chemical cleaners or abrasives for cleaning.

Important

The belt and control unit are not suitable for machine washing.

Cleaning the belt

Once you have removed the control unit from the belt, you can clean the belt using a soft cloth. To do so, use a mild soap solution or liquid detergent. Do not use bleaching agents.

8. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life.

Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

9. Troubleshooting

Problem/question	Possible cause/remedy
The unit does not switch on when you press the ON/OFF button.	<ol style="list-style-type: none"> 1. Ensure that the batteries are inserted correctly and have contact. The on/off button must be pressed and held down to switch the unit on. 2. If necessary, replace the batteries. 3. Please contact customer service.
No stimulation can be felt.	<ol style="list-style-type: none"> 1. Increase the pulse intensity incrementally. 2. Press the on/off button to restart the program. 3. The batteries are almost empty. Replace them.
You feel an unpleasant sensation from the electrodes.	The electrodes are poorly moistened. This can lead to skin irritation, as the current is no longer distributed evenly over the entire surface. Switch off the unit. Take off the belt, moisten the belt (electrodes) and your lower back and restart the program according to the manual.
The skin in the treatment area becomes red.	Discontinue the treatment immediately and wait until the skin condition has normalised. Temporary skin reddening under the electrode is harmless and can be explained by the localised increase in blood flow. However, if the skin irritation persists and you feel itching or irritation, consult your doctor before continuing use. The cause may be an allergy to the surface of the electrode.

10. Technical specifications

Belt size for waist sizes from	Approx. 75 to 140 cm
Weight (belt and unit)	Approx. 415 g with batteries
Electrode size	Approx. 81 x 65 mm
Parameters (500 Ohm load)	Output voltage: max. 50 Vpp / 5.5 Vrms Output current: max. 100 mApp / 11 mArms Output frequency: 2 – 110 Hz
Pulse duration	60 – 220 µs per phase
Wave form	Symmetric, biphasic rectangular pulses
Power supply	4.5 V (3 x 1.5 V AAA, type LR 03)
Operating conditions	0 °C to 40 °C, 20 to 65% relative humidity
Storage	0 °C to 55 °C, 10 to 90% relative humidity

The serial number is located on the device or in the battery compartment.

We reserve the right to make any technical alterations that are necessary in order to improve and develop the product further. If the unit is used outside of the specifications, proper function is no longer guaranteed!

This unit complies with European standards EN60601-1, EN60601-1-2 (In accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3 and IEC 61000-4-8) and is subject to special precautions regarding electromagnetic compatibility. Note that portable and mobile RF communication equipment can affect this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.

The unit conforms to the requirements of the European Directive for Medical Products 93/42/EC, the MPG (German Medical Devices Act). For this device, a functional test and instruction in accordance with Section 5 of the Medical Devices Operator Ordinance (MPBetreibV) is not required. It is also not necessary to carry out safety checks in accordance with Section 6 of the Medical Devices Operator Ordinance (MPBetreibV).

11. Notes on electromagnetic compatibility

Warning

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts, including all cables included in delivery. Failure to comply with the above can impair the performance of the device.
- Failure to comply with the above can impair the performance of the device.

12. Warranty/service

In case of a claim under the warranty please contact your local dealer or the local representation which is mentioned in the list "service international".

In case of returning the unit please add a copy of your receipt and a short report of the defect.

The following warranty terms shall apply:

1. The warranty period for BEURER products is either 5 years or- if longer- the country specific warranty period from date of Purchase.
In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Repairs (complete unit or parts of the unit) do not extend the warranty period.
3. The warranty shall not be valid for damages because of
 - a. improper treatment, e.g. nonobservance of the user instructions.
 - b. repairs or tampering by the customer or unauthorised third parties.
 - c. transport from the manufacturer to the consumer or during transport to the service centre.
 - d. The warranty shall not be valid for accessories which are subject to normal wear and tear (cuff, batteries etc.).
4. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

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